

Product Datasheet

Human Aldehyde Dehydrogenase 1-A3 Protein (orb1475271)

Description	Recombinant Human Aldehyde Dehydrogenase 1-A3/ALDH1A3 (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5, 20% Glycerol.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 150mM NaCl, 20% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Pro512
Source	E.coli
Uniprot ID	P47895
Expiration Date	6 months from date of receipt.