

Product Datasheet

Human Ribose-Phosphate Pyrophosphokinase 2 Protein (orb1475302)

Description	Recombinant Human Ribose-Phosphate Pyrophosphokinase 2/PRPS2 (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Protein Sequence	Pro2-Leu318
Source	Human Cells
Uniprot ID	P11908
Expiration Date	6 months from date of receipt.