

Product Datasheet

Human β -Defensin 4 Protein (orb1475328)

Description	Recombinant Human β -Defensin 4/DEFB4
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 130mM NaCl, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 130mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Glu23-Pro72
Source	E.coli
Uniprot ID	Q8WTQ1
Expiration Date	6 months from date of receipt.