

Product Datasheet

Human ILDR2 Protein (orb1475338)

Description	Recombinant Human ILDR2 (C-6His)
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Preservatives	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Protein Sequence	Leu21-Glu186
Source	Human Cells
Uniprot ID	Q71H61
Expiration Date	6 months from date of receipt.