

## Product Datasheet

### Mouse CD24 Protein (orb1475410)

<b>Description</b>	Recombinant Mouse CD24 (C-Fc)
<b>Endotoxins</b>	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
<b>Preservatives</b>	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Asn27-Gly53
<b>Source</b>	Human Cells
<b>Uniprot ID</b>	<b>P24807</b>
<b>Expiration Date</b>	6 months from date of receipt.