

Product Datasheet

Mouse Siglec-10 Protein (orb1475421)

Description	Recombinant Mouse Siglec-10 (C-Fc)
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Preservatives	Lyophilized from a 0.2 μm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Met19-Lys543
Source	Human Cells
Uniprot ID	Q80ZE3
Expiration Date	6 months from date of receipt.