

Product Datasheet

Human DIPP3-alpha Protein (orb1475465)

Description	Recombinant Human DIPP3-alpha/NUDT10(K90E)(C-His)
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	K2-P164(K90E)
Source	E. coli
Uniprot ID	Q8NFP7
Expiration Date	6 months from date of receipt.