

## Product Datasheet

### Human Translin-associated protein X Protein (orb1475470)

|                         |                                                                                                                                                                                                              |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Recombinant Human Translin-associated protein X/TSNAX                                                                                                                                                        |
| <b>Preservatives</b>    | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4                                                                                                                         |
| <b>Form/Appearance</b>  | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4                                                                                                                         |
| <b>Storage</b>          | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| <b>Tag</b>              | tag free                                                                                                                                                                                                     |
| <b>Note</b>             | For research use only                                                                                                                                                                                        |
| <b>Protein Sequence</b> | M1-S290                                                                                                                                                                                                      |
| <b>Source</b>           | E. coli                                                                                                                                                                                                      |
| <b>Uniprot ID</b>       | <b>Q99598</b>                                                                                                                                                                                                |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                                                                                                                               |