

Product Datasheet

Human Syntaxin-4 Protein (orb1475538)

Description	Recombinant Human Syntaxin-4/STX4(C-His)
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	R2-K272
Source	E. coli
Uniprot ID	Q12846
Expiration Date	6 months from date of receipt.