

Product Datasheet

Human Leucine-rich repeat-containing protein 59 Protein (orb1475574)

Description	Recombinant Human Leucine-rich repeat-containing protein 59/LRRC59(C-His)
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	T2-S244
Source	E. coli
Uniprot ID	Q96AG4
Expiration Date	6 months from date of receipt.