


Product Datasheet

c-Jun (Phospho-Tyr170) Antibody (orb14757)

Description	c-Jun (Phospho-Tyr170) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05412
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...