

Product Datasheet

Human Nuclear apoptosis-inducing factor 1 Protein (orb1475745)

Description	Recombinant Human Nuclear apoptosis-inducing factor 1/NAIF1(C-His)
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	A2-L282
Source	E. coli
Uniprot ID	Q69YI7
Expiration Date	6 months from date of receipt.