

Product Datasheet

Human Uracil phosphoribosyltransferase homolog Protein (orb1475762)

Description	Recombinant Human Uracil phosphoribosyltransferase homolog/UPRT(C-His)
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	A2-A190
Source	E. coli
Uniprot ID	Q96BW1
Expiration Date	6 months from date of receipt.