

Product Datasheet

Human Ragulator complex protein LAMTOR3 Protein (orb1476138)

Description	Recombinant Human Ragulator complex protein LAMTOR3(C-His)
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 500mM NaCl, 15%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 500mM NaCl, 15%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	M1-S124
Source	E. coli
Uniprot ID	Q9UHA4
Expiration Date	6 months from date of receipt.