



Product Datasheet

ATF2 (Phospho-Ser62 or 44) Antibody (orb14763)

Description	ATF2 (Phospho-Ser62 or 44) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P15336
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of ATF2(Phospho-Se...