

Product Datasheet

Human TFB1M Protein (orb1476306)

Description	Recombinant Human TFB1M Protein(N-His-SUMO)
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	N-His-SUMO
Note	For research use only
Protein Sequence	Y58-S194
Source	E. coli
Uniprot ID	Q8WVM0
Expiration Date	6 months from date of receipt.