

Product Datasheet

Human Polycomb group protein ASXL1 (orb1476341)

Description	Recombinant Human Polycomb group protein ASXL1
Preservatives	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	H61-R246
Source	E. coli
Uniprot ID	Q8IXJ9
Expiration Date	6 months from date of receipt.