

Product Datasheet

Human FcRn & B2M Heterodimer Protein (orb1476378)

Description	Recombinant Human FcRn & B2M Heterodimer
Endotoxins	Less than 0.1 EU/ug
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, 0.05%Tween 80, pH 7.4
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, 0.05%Tween 80, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Protein Sequence	Ala24-Ser297&Ile21-Met119
Source	Human Cells
Uniprot ID	P55899, P61769
Expiration Date	6 months from date of receipt.