

Product Datasheet

ERp57 (GRP58) antibody (orb147640)

Description

Mouse monoclonal to Erp57 (PerCP). ERp57,

(Grp58), Hormon...

Species/Host

Mouse

Reactivity

Bovine, Canine, Guinea pig, Hamster,
Human, Monkey, Mouse, Porcine, Rabbit,
Rat

Conjugation

PerCP

Tested

ICC, IHC, IP, WB

Applications
Immunogen

Human recombinant ERp57 (Grp58)

Target

ERp57

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium
azide *Storage buffer may change when
conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at
4°C

Note

For research use only

Application notes

0.5 µg/ml of SMC-168 was sufficient for
detection of ERp57 in 10 µg of heat shock
Heal Lysate by colorimetric immunoblot
analysis using Goat anti-mouse IgG:HRP as
the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

57kDa

Uniprot ID

[P30101](#)

NCBI

[NP_005304.3](#)

Entrez

[2923](#)

Dilution Range

WB (1:2000), IHC (1:100), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.