

Product Datasheet

Human IGFBP-7 Protein (orb1476524)

Description	Recombinant Human IGFBP-7 (C-His)
Endotoxins	Less than 0.1 EU/ug
Preservatives	Lyophilized from a 0.2 um filtered solution of PBS, 0.05%Tween 80, pH 7.4
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, 0.05%Tween 80, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Protein Sequence	Ser27-Leu282
Source	Human 293 Cells
Uniprot ID	Q16270
Expiration Date	6 months from date of receipt.