

Product Datasheet

AHA1 antibody [Out of stock] (orb147660)

Description

Rat monoclonal to Aha1 (ATTO-390). Aha1

family, and is thought ...

Species/Host

Rat

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-390

Tested

ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Recombinant Full Length Mouse Aha1 Protein

Target

AHA1

Preservatives

0.02M potassium phosphate, 0.15M sodium chloride, pH 7.2 *Storage buffer may change when conjugated

Concentration

1.2 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-172 was sufficient for detection of Aha1 in 10 µg of rat tissue lysate by colorimetric immunoblot analysis using Goat anti-rat IgG:HRP as the secondary antibody.

Isotype

IgG2a Kappa

Clonality

Monoclonal

MW

38kDa

Uniprot ID
[Q8BK64](#)
NCBI
[NP_666148.1](#)
Entrez
[217737](#)
Dilution Range

WB (1:1000), IHC (1:100), ICC/IF (1:1000), IP (1:1000)

Expiration Date

12 months from date of receipt.