

Product Datasheet

Human TREM-1 Protein (orb1476627)

Description	Recombinant Human TREM-1 (C-His)
Endotoxins	Less than 0.1 EU/ug
Preservatives	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Protein Sequence	Ala21-Arg200
Source	Human 293 Cells
Uniprot ID	Q9NP99
Expiration Date	6 months from date of receipt.