

Product Datasheet

AHA1 antibody (orb147674)

Description

Rat monoclonal to Aha1 (PerCP). Aha1 is a

is thought to ...

Species/Host

Rat

Reactivity

Human, Mouse, Rat

Conjugation

PerCP

Tested

ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Recombinant Full Length Mouse Aha1 Protein

Target

AHA1

Preservatives

0.02M potassium phosphate, 0.15M sodium chloride, pH 7.2 *Storage buffer may change when conjugated

Concentration

1.2 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-172 was sufficient for detection of Aha1 in 10 µg of rat tissue lysate by colorimetric immunoblot analysis using Goat anti-rat IgG:HRP as the secondary antibody.

Isotype

IgG2a Kappa

Clonality

Monoclonal

MW

38kDa

Uniprot ID

[Q8BK64](#)

NCBI

[NP_666148.1](#)

Entrez

[217737](#)

Dilution Range

WB (1:1000), IHC (1:100), ICC/IF (1:1000), IP (1:1000)

Expiration Date

12 months from date of receipt.