


Product Datasheet

STAT1 (phospho-Tyr701) antibody (orb14776)

Description	STAT1 (phospho-Tyr701) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human STAT1 around the phosphorylation site of Tyr701 (TG(p-Y)IK)
Target	phospho-STAT1 (Tyr701)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82 kDa
Uniprot ID	P42224
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



From [Cancer Medicine](#) (2016.6):
Publ...