

## Product Datasheet

### Acrosin antibody (orb1478)

## Description

Acrosin antibody

### Species/Host

Rabbit

### Reactivity

Human, Rat

### Conjugation

Unconjugated

### Tested Applications

ICC, IF, IHC-P, WB

### Immunogen

KLH conjugated synthetic peptide derived from human Acrosin. Please contact us for the exact immunogen sequence. The peptide is available as orb393925.

### Target

Acrosin

### Form/Appearance

10 mM PBS, 0.02% sodium azide

### Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

### Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

### Note

For research use only

### Isotype

IgG

### Clonality

Polyclonal

### Purity

Polyclonal antibodies are purified by peptide affinity chromatography

### MW

45 kDa

### Uniprot ID

[P10323](#), [P29293](#)

### NCBI

[001088](#), [001097](#)

### Entrez


[49](#)

### Dilution Range


IF/ICC: 1:50-300, IHC-P: 1:50-300, WB: 1:200-1:800

### Expiration Date


12 months from date of receipt.



WB analysis of rat ovary (lane 1), mouse...



IHC-P image of human breast cancer tissue...



Immunofluorescence analysis of human breast cancer tissue...