


Product Datasheet

ATF4 (Phospho-Ser245) Antibody (orb14784)

Description	ATF4 (Phospho-Ser245) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P18848
Expiration Date	12 months from date of receipt.



P-Peptide

Immunohistochemical analysis of paraffin...