

Product Datasheet

CD45 antibody (orb1478764)

Description	Mouse monoclonal antibody to CD45
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P, IP
Immunogen	Human T lymphocytes
Preservatives	PBS with 0.1% sodium azide
Form/Appearance	Purified
Concentration	IgG 1 mg/ml
Storage	aliquots -20°C. Thawed 2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.