

Product Datasheet

CD29 antibody (orb1478769)

Description	Mouse monoclonal antibody to CD29
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, IP
Immunogen	Occular melanoma cell line V+B2
Preservatives	PBS with 0.1% sodium azide
Form/Appearance	Purified
Concentration	IgG 1 mg/ml
Storage	aliquots -20°C. Thawed 2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.