

Product Datasheet

TIA antibody (orb1478801)

Description	Mouse monoclonal antibody to TIA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Human bone marrow malignant cells from a non-B, non-T acute leukemia
Preservatives	Liquid purified Ascites, purified with Protein G Chromatography with 15mM Sodium Azide
Form/Appearance	Purified
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC-P (1:100-1:200)
Expiration Date	12 months from date of receipt.