

Product Datasheet

P16 antibody (orb1478809)

Description	Mouse monoclonal antibody to P16
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Unknown or proprietary to Biorbyt and/or its suppliers
Preservatives	Tissue culture supernatant containing 15mM sodium azide
Form/Appearance	Concentrate
Storage	2-8°C
Note	For research use only
Isotype	IgGk
Clonality	Monoclonal
Dilution Range	IHC-P (1:100-1:400)
Expiration Date	12 months from date of receipt.