

Product Datasheet

NOD2 (C) antibody (orb1478845)

Description	Rabbit polyclonal antibody to NOD2 (C)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NOD2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	115
Uniprot ID	Q9HC29
Expiration Date	12 months from date of receipt.