

Product Datasheet

Zika E (M) antibody (orb1478849)

Description	Rabbit polyclonal antibody to Zika E (M)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of Zika virus (strain Mr 766) envelope protein E.
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	54
Uniprot ID	Q32ZE1
Expiration Date	12 months from date of receipt.