

Product Datasheet

RAD51 (C) antibody (orb1478899)

Description	Rabbit polyclonal antibody to RAD51 (C)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of arabidopsis thaliana RAD51
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	37
Uniprot ID	P94102
Expiration Date	12 months from date of receipt.