

Product Datasheet

CLV1 (N) antibody (orb1478959)

Description	Rabbit polyclonal antibody to CLV1 (N)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of arabidopsis thaliana CLV1 (At1g75820)
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	108
Uniprot ID	Q9SYQ8
Expiration Date	12 months from date of receipt.