

Product Datasheet

RBCL (C) antibody (orb1478978)

Description	Rabbit polyclonal antibody to RBCL (C)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of arabidopsis thaliana RBCL (AtCg00490)
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	53
Uniprot ID	003042
Expiration Date	12 months from date of receipt.