

Product Datasheet

TOC1 (C) antibody (orb1478998)

Description	Rabbit polyclonal antibody to TOC1 (C)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of arabidopsis thaliana TOC1 (At5g61380)
Preservatives	0.01% NaN ₄
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	69
Uniprot ID	Q9LKL2
Expiration Date	12 months from date of receipt.