

Product Datasheet

Piccolo antibody [Out of stock] (orb147918)

Description

Mouse monoclonal to Piccolo (ATTO-594).

	large protein whi...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-594
Tested Applications	ELISA, ICC, WB
Immunogen	Full length recombinant protein (GenBank ID: BC001304.1)
Target	Piccolo
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-188 was sufficient for detection of Piccolo in 20 µg of rat brain tissue lysate by colorimetric immunoblot analysis using goat IgG:HRP as the secondary antibody.
Isotype	IgG2b Kappa
Clonality	Monoclonal
MW	550kDa
Uniprot ID	Q9QYX7
NCBI	NP_001104266.1
Entrez	26875
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.