

Product Datasheet

TAT-HuR-HNS3 (orb1479252)

Description	HuR-PARP1 interaction blocker; Cell penetrating peptides.
Form/Appearance	Freeze dried solid
Storage	Store dry, frozen and in the dark
Note	For research use only
Protein Sequence	YGRKKRRQRRSPMGVDHMSGVNVPGNASSG
Formula	$C_{151}H_{257}N_{59}O_{46}S_2$
Purity	> 95% by hplc
MW	3696.9 Da
Solubility (25°C)	Soluble in water
Expiration Date	12 months from date of receipt.