

Product Datasheet

Human ULBP3 ELISA Kit (orb1479262)

Description	Human UL16-binding protein 3(ULBP3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: tissue homogenates, cell lysates and other biological fluids; sandwich
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9BZM4
Sensitivity	0.108ng/mL
Expiration Date	Please enquire.