

Product Datasheet

Human ULBP1 ELISA Kit (orb1479263)

Description	Human UL16-binding protein 1(ULBP1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: tissue homogenates, cell lysates and other biological fluids; sandwich
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9BZM6
Sensitivity	0.058ng/mL
Expiration Date	Please enquire.