

Product Datasheet

Mouse OXTR ELISA Kit (orb1479282)

Description	Mouse Oxytocin Receptor (OXTR) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: tissue homogenates, cell lysates and other biological fluids; sandwich
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P97926
Sensitivity	0.048ng/mL
Entrez	18430
Expiration Date	Please enquire.