

Product Datasheet

Mouse DEFb3 ELISA Kit (orb1479304)

Description	Mouse Defensin Beta 3 (DEFb3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: serum, plasma and other biological fluids; sandwich: Experimental Method
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.052ng/mL
Expiration Date	Please enquire.