

Product Datasheet

Mouse CTRC ELISA Kit (orb1479313)

Description	Mouse Chymotrypsin-C (CTRC) ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: serum, plasma and other biological fluids; sandwich: Experimental Method
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q3SYP2
Sensitivity	0.34ng/mL
Entrez	76701
Expiration Date	Please enquire.