

Product Datasheet

Human CTRC ELISA Kit (orb1479314)

Description	Human Chymotrypsin-C (CTRC) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: serum, plasma, tissue homogenates and other biological fluids; sandwich: Experimental Method
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q99895
Sensitivity	1.35ng/mL
Entrez	11330
Expiration Date	Please enquire.