

Product Datasheet

Human C1qB ELISA Kit (orb1479326)

Description	Human Complement Component 1, Q Subcomponent B (C1qB) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: serum, plasma and other biological fluids; sandwich
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02746
Sensitivity	0.055ng/mL
Entrez	713
Expiration Date	Please enquire.