

Product Datasheet

Mouse AGR2 ELISA Kit (orb1479328)

Description	Mouse Anterior Gradient Protein 2 (AGR2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids; sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O88312
Sensitivity	0.059ng/mL
Entrez	23795
Expiration Date	Please enquire.