

Product Datasheet

Monkey aFP ELISA Kit [Out of stock] (orb1479329)

Description	Monkey Alpha-Fetoprotein (aFP) ELISA Kit
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Note	For research use only
Application notes	Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids; sandwich: Experimental Method
Sensitivity	1.2ng/mL
Expiration Date	12 months from date of receipt.