

Product Datasheet

Mouse XCL1 ELISA Kit (orb1479332)

Description	Mouse Chemokine C-Motif Ligand 1 (XCL1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: tissue homogenates, cell lysates and other biological fluids; sandwich
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P47993
Sensitivity	0.055ng/mL
Entrez	16963
Expiration Date	Please enquire.