

## Product Datasheet

### Mouse OXTR ELISA Kit (orb1479356)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Oxytocin Receptor (OXTR) ELISA Kit                                                 |
| <b>Reactivity</b>        | Mouse                                                                                    |
| <b>Range</b>             | 0.156-10ng/mL                                                                            |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | Sample Type: tissue homogenates, cell lysates and other biological fluids; sandwich      |
| <b>Sample Types</b>      | tissue homogenates, cell lysates and other biological fluids.                            |
| <b>Uniprot ID</b>        | <b>P97926</b>                                                                            |
| <b>Sensitivity</b>       | 0.048ng/mL                                                                               |
| <b>Entrez</b>            | <b>18430</b>                                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                                          |