

Product Datasheet

Human METRNL ELISA Kit (orb1479362)

Description	Human Meteorin-Like Protein (METRNL) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: serum, plasma, tissue homogenates and other biological fluids; sandwich
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q641Q3
Sensitivity	12.8pg/mL
Entrez	284207
Expiration Date	Please enquire.